

Pro Ligation-Free Cloning Kit: Quick-Start Protocol

To download the full manual, visit Cat. No. E086 or E087 on www.abmgood.com.



Design and Amplify Overlapping PCR Products

- 1. Design primers to amplify insert fragments with 15-25 bp of overlapping homology to each other and to the linearized vector.
- 2. PCR amplify the fragments and visualize on an agarose gel.

Critical Step: For more details on primer design, see page 4 of the manual.

Purification Steps (Recommended)

- 3. Column purify your products, or gel extract if the PCR amplification contained non-specific bands.
- Alternatively, add 0.5 µl of Assembly Enhancer and 1.0 µl Cloning Optimizer to 50 µl of PCR product. Incubate at 37°C for 5 minutes, then at 80°C for 20 minutes.

Note: Gel extraction will help remove any residual PCR template plasmid, reducing background colony growth.

Assembly Reaction

- 4. Quantify the insert and vector DNA by a UV spectrophotometer or by comparing the target band against a known molecular weight marker run on the same gel.

- 5. Set up the following reaction on ice:

Reagent	Cloning Reaction	Example Reaction
2X Pro Ligation-Free MasterMix	10 µl	10 µl
Linearized Vector	100 ng	100 ng (5 kb vector)
Each Insert	3:1 insert-vector molar ratio	60 ng (1 kb insert)
Nuclease Free H ₂ O	Up to 20 µl	Up to 20 µl

To calculate amount of insert to add:

$$\text{insert (ng)} = 3 \times \frac{\text{insert (bp)}}{\text{vector (bp)}} \times \text{vector (ng)}$$

- 6. For 1 insert, incubate at 50°C for 15 minutes. For multiple inserts or difficult assemblies, incubate at 50°C for 1 hour.

- 7. Transform immediately, or store samples at -20°C until transformation.

ProClone™ Chemical Transformation

- 8. Add reaction mix to 60 µl of ProClone™ competent cells (Cat. No. E003).
- 9. Incubate on ice for 30 minutes.
- 10. Heat-shock for 45 seconds at 42°C, then chill on ice for 2 minutes.
- 11. Add 150 µl of LB Medium without antibiotics.
- 12. Shake at 37°C for 1 hour for recovery.
- 13. Spread cells onto pre-warmed LB agar plates with an appropriate antibiotic.
- 14. Incubate plates overnight at 37°C.

Troubleshooting: If no colonies are present, check: primer design, purity and concentration of DNA inserts, vector-insert ratio, antibiotic used, and transformation efficiency.

Screening and Analysis

- 15. Screen by colony PCR, or grow up colonies overnight in LB medium containing the appropriate antibiotic.
- 16. Isolate plasmid DNA. We recommend using our Plasmid Mini-Prep Kit (D504).
- 17. Confirm correct assembly by restriction analysis or sequencing.

Troubleshooting: If too many colonies are present or clones do not pass screening, consider the following possibilities: vector was not completely digested, presence of residual PCR template plasmid, antibiotics are old, or amplification was not specific.